

## CT26-hTROP2

## **Strain Information**

- SMOC番号 NM-S24A-TM38
- Cell Line CT26-hTROP2

維持形態 Validation of tumorigenic capacity completed

説明 TROP2 gene humanized CT26 cell line

\*Literature published using this strain should indicate: CT26-hTROP2 cell line (Cat. NO. NM-S24A-TM38) was purchased from Shanghai Model Organisms Center, Inc..

## 表現型デロタ




Fig.1 hTROP2 expression on in vitro cultured cells of CT26-hTROP2.






Fig.2 In vivo validation of CT26-hTROP2 in BALB/c mice.BALB/c mice were inoculated with CT26hTROP2 subcutaneously at 3×106 cells/mouse inoculum for tumor development. The day of vaccination for mice is defined as Day 0, and tumor measurement and weight measurement begin on the 7th day after vaccination, twice a week.




Fig.3 hTROP2 expression in tumors after in vivo growth of CT26-hTROP2 in wild type BALB/c mice.When the MTV reaches 500-800 mm<sub>3</sub>, take one tube from each group and send it to FACS for testing (Markers: L/D, mCD45, hTROP2). A, the MTV reaches 800 mm<sub>3</sub>. B, the MTV reaches 1500 mm<sub>3</sub>.